Afname Hemocultuur.

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1. BD BACTEC vial & skin preparation

- Label bottle with collection date, time, and patient information (do not cover bottle barcode label)

- Disinfect the skin (>60 seconds). Remove plastic flip-off cap from the vial and disinfect the rubber septum with a single alcohol swab and allow to dry.

- Mark bottle label corresponding to the recommended fill level.

2a. Collecting Blood Sample Using BD Vacutainer® Safety-Lok™ Blood Collection Set

- Using **BD Vacutainer® Safety-Lok Blood Collection Set** and holder perform venipuncture by holding wings as shown. **DO NOT** hold by grasping the safety shield.

- Push and hold Vacutainer holder over top of vial to puncture septum.

- First inoculate blood culture vials before drawing blood in regular Vacutainer blood collection tubes. Fill Plus Aerobic bottle first, then Lytic Anaerobic bottle. Hold the bottle upright. Collect blood to indicated fill level on vial. Recommended blood volume is 8-10 ml. Accepted volume is 3-10 ml. For pediatric vials recommended blood volume is 0,5-3 ml. For Myco F/Lytic recommended blood volume is 5 ml.

- If required, additional regular Vacutainer tubes may be drawn at this time by inserting them into the Vacutainer holder.

2b. Collecting Blood Sample Using BD Vacutainer® Catheter Blood Collection Set

- Connect the Catheter Blood Collection set to the Single Use Holder.

- Insert the male luer adapter into the catheter hub. Secure the Connection with a 1/4 turn clockwise.

- Collect the blood, respecting the recommended order as explained in section 2a.

- Once sampling is complete, stabilize the catheter and gently pull back the catheter blood collection set applying a 1/4 reverse clock turn.

3. Remove needle from venipuncture site

- When the final tube is filled, withdraw the needle by grasping the translucent yellow safety shield grip area. **DO NOT** withdraw by holding the safety shield.

- With opposite hand grasp tubing between thumb and index finger. To activate the safety shield, push the yellow shield forward until it is locked in place. Dispose the nearest sharps container without disassembling.

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